  
BRASTER<sup>®</sup>

# About Braster

**Braster S.A.** is an innovative Polish telemedical company listed on the Warsaw Stock Exchange. It was founded in 2008 by a group of scientists who have developed a state-of-the-art technology called contact thermography.


In October 2016 Braster S.A. launched **Braster System** - innovative medical device for in-home breast self-examination.

In August 2018 we launched **Braster Pro** – innovative medical device for professional use by Healthcare professionals.

# Production plant


# Quality control

# Cancer growth

There are two major biological processes that influence cancer growth:

- proliferation
- angiogenesis


Angiogenesis is the process of forming new vessels, and it plays an essential role in the development of breast cancer and its metastasis


# Cancer growth

These two processes result in an increase in surrounding tissue metabolism, and subsequently a higher temperature in the place where the cancer is growing

The increase in temperature starts in the very early stages of cancer development!



BRASTER<sup>®</sup>

# Method and Product


**Thermography**- imaging, detection and registration of the temperature in the examined body.

**Remote-sensing thermography** - the image is obtained without the need of contact of the device and the Surface of examined body; the technique is based on heat transfer through radiation [e.g. infrared (IR) thermographic cameras, pyrometers];


**Contact thermography** - the image is obtained through the contact of the device with examined surface based on heat transfer by thermal conduction [e.g. liquid-crystal forehead thermometers, Braster device].


# Unique liquid-crystal technology


Liquid crystals are chemical compounds that exhibit the properties of liquids and those of crystalline solids. BRASTER S.A. has an innovative, proprietary technology for producing mixtures of liquid crystals and a Continuous Liquid Crystal Film (CLCF) technology for applying liquid-crystal emulsion on the film.


# Thermographic matrices

BRASTER PRO is delivered with a set of 3 thermographic matrices which use liquid crystals to present breast surface heat distribution in the form of a colorful map



Each matrix in a set is calibrated to a different temperature range to compensate for differences in body temperature across patients:


# Image acquisition device

Heat map presented on the surface of liquid crystal matrix during examination is illuminated with LEDs and captured with visible light spectrum camera

# Main components

Housing


Hand Grip / Examination button

Digital camera

Li-Ion rechargeable battery


Main PCB

LED strip

Light chamber

Light shield, diffuser

# Braster Pro breast examination

## Mobile app

guides through the examination  
and informs about the result


Device  
Braster Pro

Images are sent for  
interpretation to  
telemedical center

Certified medical experts  
control of the interpretation  
process


Examination result

Telemedical center  
where the automatic  
interpretation of the  
breast thermographic  
images is performed  
(Braster AI)

# System for automatic interpretation Braster AI

New generation algorithms - **DEEP LEARNING NETS**


# Braster AI Algorithms

Thermal images of two breasts are analysed and compared in three main parameters:

**Thermal Asymmetry** – a thermal difference between hottest areas in both breasts

**Structural Asymmetry** – difference in number of the thermal structures between both breasts

**Area Asymmetry**- a significant disproportion between the areas of main structures of both breasts

Details of Braster® AI analysis

	Value:	Norm:
Thermal asymmetry parameter:	0.00	from -2 to 2
Structure asymmetry parameter:	0.00	from -2 to 2
Main thermal structure area ratio:	-3.12	from -4 to 4
Main difference ratio of surface thermal structures areas:	-42216.50	from -20,000 to 20,000


# Braster Pro – clinical validation

**1. BRA/03/2013 (ThermaCRAC)** – „Multicentre, observational study comparing the effectiveness of the Tester BRASTER™ device in diagnosing and differentiating breast pathology in women to standard diagnostic methods; n = 736 (proof of concept)

**2. BRA 03/2014 (ThermaRAK)** - multicentre observational study conducted to collect data from thermographic examination and data from imaging and histopathological examinations necessary to prepare the atlas of thermal pathologies. n = 360;

**3. BRA/11/2014 (ThermaALG)** – Multicentre observational study assessing the diagnostic effectiveness and clinical usefulness of the new version of the interpretation algorithm of the thermographic examination in the diagnosis of breast pathology in women; n=274 (registration study)

**4. INNOMED\_BRASTER\_2016\_01**– Multicentre observational study assessing diagnostic effectiveness of contact thermography in comparison with ultrasound examination, mammography and breast biopsy; n=3000 (PMS)

# Observational study ThermaAlg

<https://clinicaltrials.gov/ct2/show/NCT03858738>

A multicenter study conducted to assess the diagnostic effectiveness of the interpretation algorithm of thermographic images in detecting breast cancer.

**Methodology:** the study included 3 groups of patients, two symptomatic groups with USG or MMG score BIRADS 4/5 divided according to age into patients (50- and 50+) and a control group with the USG score BIRADS 1/2 (n=274)

**The primary goal** of the study was to determine the effectiveness of diagnostic contact thermography using the manual algorithm of assessment the thermographic images compared to breast ultrasound, mammography and histopathological examination of the lesion.

**The secondary goal** was the validation of automatic algorithms using artificial intelligence in the assessment of thermographic images.

## Patients under 50 years of age:

- ❖ Sensitivity 81,5 % (95% CI 64,1; 92,6)
- ❖ Specificity 87% (95%CI: 79,7; 92,4)
- ❖ PPV 71,0% (95%CI: 53,7; 85,8)
- ❖ NPV 93,6% (95%CI: 85,6; 97,8)

# ThermaALG (clinical trial) conclusions

## Braster + US

💡 **positive** result in contact thermography **increases** the likelihood of diagnosis of breast cancer more than **twofold**

💡 **negative** result in contact thermography **decreases** the likelihood of diagnosis of breast cancer more than **threefold**

# Observational study INNOMED

<https://clinicaltrials.gov/ct2/show/NCT03858738>

Prospective multicentre (**24 sites**) observational study conducted by Collegium Medicum Jagiellonian University in Krakow.

**Methodology:** 3000 patients; 3 groups : A: women US BIRADS 4-5 < 50 C: women US BIRADS 4-5 ≥ 50 B: control group of women BIRADS 1-2; 18-49; ≥ 50

## The primary goal :

Diagnostic effectiveness of contact thermography in comparison with ultrasound examination, mammography and breast biopsy

## The secondary goal :

Validation of algorithms for automatic interpretation of thermographic images (deep learning algorithm)

## Results Q4 2019

- The average sensitivity of contact thermography is 60.7-61.8%
- The specificity of contact thermography is between 62.7 – 85.4%
- The results obtained indicate the limited use of contact thermography as an independent screening method

# INNOMED Study Results


- The average **sensitivity** of contact thermography is 60.7-61.8%
- The **specificity** of contact thermography is between 62.7 – 85.4%

PPV in the group of women <50

BIRADS	PPV-usg	CNT-usg	PPV-Braster+	CNT-Braster+	PPV-Braster-	CNT-Braster-
4a	-	0	-	0	-	0
4b	8.5%	71	19.2%	26	2.2%	45
4c	75.0%	8	100.0%	5	33.3%	3
4a-c	15.2%	79	32.3%	31	4.2%	48
5	100.0%	14	100.0%	10	100.0%	4

ROC Curve (sensitivity-specificity)

Braster + US

BIRADS 4 (from ultrasonography) with a positive result in **thermography increases the positive predictive value more than twofold**, while a negative result of this test **significantly reduces** this value.

This confirms the **use of contact thermography as a complementary method** to breast ultrasound (similar results were obtained in the BRA11 / 2014 ThermaAlg study)

# Braster Pro - intended use & contraindication

- is intended for breast examination of women over 18 by HCP
- is an adjunct to recognized modalities such as ultrasound and mammography
- detects thermal asymmetry in women's breasts which can be correlated with breast pathology
- proved efficient regardless of the size of the breast, its density or aesthetic implants

## **Contraindications** for Braster Pro examination:

- **Patient who are undergoing or have completed anti-cancer therapy of breast cancer**
- Temporary contraindications:
  - pregnancy or breastfeeding (up to three months after weaning);
  - general infection, with a body temperature of or in excess of 38°C;
  - breast infection with pain, skin redness and bruises (when said symptoms are present);
  - surgical procedure in the breast area with benign lesion diagnosis:
    - ❖ fine-needle biopsy (FNB) – up to four weeks after the procedure,
    - ❖ core-needle biopsy (CNB) or Mammotome's breast biopsy – up to 6 months after the procedure,
    - ❖ breast tumor resection – up to 12 months after the procedure;
- **aesthetic implant placement**, filler injections (e.g. hyaluronic acid) and lipotransfer – **up to 12 months after the procedure.**

# Algorithm for Braster Pro

**18-49**

negative

routine check in  
planned time

positive

diagnostic  
examination

\*Polish medical societies recommend routinely performing a breast ultrasound examination at an annual time interval

**50 +**


MMG  
screening

After 2 years next  
screening MMG

# Algorithm for Braster Pro


**+** positive  
**-** negative



\*BI-RADS 4a – follow up in 6 months according to BI-RADS procedure

CBE – Clinical breast examination (PALPATION)  
US- Ultrasound

# Braster Pro- who can perform an examination

  
BRASTER<sup>®</sup>  
**Cases**


## CASE 1

**Age:** 49 years


**Breast composition:** fatty-glandular

**Focal change:** palpable


**Hist-pat:** invasive ductal carcinoma (*carcinoma ductale invasivum*)

In the left breast, an irregular focal hyperthermia which is corresponding to a verified cancer.

Ultrasound: In the left breast, at the perimeter of the 2 o'clock position, an irregular hypoechoic mass (measuring 16 x 12 mm) is visible; BI-RADS 4b.

MMG MLO: No suspicious focal changes or clusters of microcalcifications are visible; BI-RADS 1.


## CASE 2

**Age:** 39 years


**Breast composition:** dense glandular

**Focal change:** non-palpable


**Hist-pat:** invasive ductal carcinoma (*carcinoma ductale invasivum*)

A focal hyperthermia was observed in the lower outer quadrant of the left breast, in the location previously visualized through ultrasound examination.

Ultrasound examination of the left breast revealed a 17 x 24 mm oval, hypoechoic lesion, with indistinct margins in the lower outer quadrant, at 4 o'clock axis, BI-RADS 4b

Due to the type of breast tissue (i.e., dense glandular tissue according to Wolfe's classification), mammographic examination had reduced sensitivity. The mammogram did not detect any suspicious areas or clusters of calcifications. BI-RADS 0


### CASE 3

**Age:** 40 years


**Breast composition:** dense glandular

**Focal change:** non-palpable


**Hist-pat:** invasive ductal carcinoma(*carcinoma ductale invasivum*)

Thermography showed irregular hyperthermy localised in the centre of right breast, on the border of the upper quadrants.

On ultrasonography several simple cyst were noted in both breast. In the right breast at the 12 o'clock axis irregular hypoechoic lesion 33x18mm was identified, suspicious of malignancy, BI-RADS 4b

Oval, well circumscribed foci up to 10mm in diameter were seen, no suspicious mass or clusters of calcifications was found, qualified patient for an additional ultrasound evaluation, BI-RADS 0

BRASTER<sup>®</sup>

# Opinions & Publications

# Opinions

Medical center Maichin Dom

2, Zdrave str.

1463 Sofia, Bulgaria

To whom it may concern,

This is to confirm that Braster Pro system, manufactured by Braster SA in Poland, has been supplied to our Clinic Maichin Dom in Sofia.

Braster Pro is an innovative solution for breast examinations performed by healthcare professionals in clinical setting. The system consists of a high-tech device which uses contact thermography, user-friendly mobile application, automatic interpretation algorithms detecting potentially dangerous changes in breasts and telemedical centre with experienced medical experts. The examination result is available within 48 hours and is sent directly to the doctor's account. The functionality of the account allows the doctor to review all the patients' results as well as check and download thermographic images.

We want to emphasize that the patients are willing to undergo Braster examination because it is painless, safe and takes only 20 minutes. The system which includes mobile app for smartphones and tablets connects wirelessly to the device, providing intuitive instructions throughout the entire examination.

The results are credible and in case when needed they can be confirmed with other modalities.

We decided to use Braster Pro system as we strongly believe that combining various types of examination increases the chance of detecting alarming changes at an early stage of the disease, which may allow for less invasive treatment and enhance the probability of a complete recovery.

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## OPINION REGARDING THE USABILITY OF BRASTER PRO SYSTEM AS A COMPLEMENTARY METHOD FOR BREAST CANCER DIAGNOSTICS IN THE CLINICAL SETTINGS

### Introduction

Cancer is one of the most common causes of death in Poland and in the world. According to the National Cancer Registry data, in 2015, 23.6% of women died because of cancer. In 14.1% the cause of death in women was breast cancer. In 2015, the mortality rate was 14.1% in Poland and 14.1% in the world. In order to provide the best accuracy and benefit from the latest methods of treatment, it is necessary to ensure the widest possible access to effective secondary prevention.


According to the report from the Supreme Audit Office, from January 2018: "in Poland there is no comprehensive, consistent and functional system of secondary health prevention, including planning activities in this area, supervision of their implementation and evaluating the obtained effect." In 2017, doctors providing basic health services implemented limited health prevention tasks to their practice. In the breast examination, we have information about the examination or its refusal in 13.5% of patients. We can therefore assume that 86.5% of women were not a breast examination at the general practitioner's clinics.

According to a survey conducted by the Chair and Department of Preventive Health, Health Sciences, at Poznan University of Medical Sciences, the most common knowledge and preventive activities for women is from a gynecologist. The above data, as oncological-gynecologists, we face the question of how to do of early detection of breast cancer, an important element of prophylaxis, sensitivity (21-41% depending on the patient's age), is self-examination. I do not self-examine their breasts because they are not aware of how to do

## MODERN PROBLEMS REQUIRE MODERN SOLUTIONS NEW TECHNOLOGY IS SHAPING HEALTHCARE

You must have seen this statement quite a few times while scrolling through social media. The statement is quite self-explanatory, meaning that more contemporary the problem, more new-age solutions are required to counter and solve it. According to a report by PwC, in April 2019, India saw its population rise with approximately 400 million inhabitants over 25 years leading up to 2016. In the same period, the population living below the poverty line decreased by half. Here is the flipside. Due to this situation, India is facing what media outlets are calling as a "dual disease burden", which is characterized by an acceleration in communicable and non-communicable diseases, also known as lifestyle diseases. Not only in India, but we are encountering new forms of diseases, infections, and afflictions all around the world, which calls for advanced and modernized healthcare solutions.


**Vivek Tiwari**  
Founder, and CEO,  
Medikabazaar

Technology's integration with healthcare is and has been providing radical solutions to diagnose and treat patients suffering from complex afflictions. As of this moment, technologically advanced healthcare solutions are no longer in its infancy, 2020 is right around the corner and we are in an era where healthcare is taking the next big step, and this is happening because of various forms of technology. It is to be noted here that as patient diagnosis but also directed towards improving the operational efficiency of healthcare facilities. Let's look at some of the trending technologies in healthcare at the moment.

Blockchain is one technology which is revolutionizing the way hospitals collect data. Patient and transaction data is one of the most sensitive information for hospitals. Via blockchain's digital ledgers and record-keeping systems, such critical data will be secure in an age where data privacy is facing some serious questions. Virtual Reality (VR), according to reports, is poised for a market valuation of \$4 billion by 2020. This technology can help hospitals create new ways of engaging with their patients (virtual tour of facilities, help patients manage their pain), which can help them find less invasive about visiting medical establishments.

I am keeping the best for last. Artificial Intelligence (AI) is the technology which is not only being implemented in you

Passes and safety stock indicator page on VUI

WHY BRASTER<sup>®</sup> PRO ?

# Braster Pro

-  1. Combining various methods increases the efficiency of breast examination
-  2. Certified medical device – CE certificate
-  3. Effectiveness proven in observational studies
-  4. Detects non - palpable lesions (the smallest lesion detected was 3mm)
-  5. Painless, radiation free and efficient regardless of breast density and size
-  6. Handy, mobile, easy to deisinfest and keep clean
-  7. Affordable and easily available

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2019 KBS COPYRIGHT KOREAN BREAST CANCER SOCIETY. ALL RIGHTS RESERVED.  
2019 KBS COPYRIGHT KOREAN BREAST CANCER SOCIETY. ALL RIGHTS RESERVED.  
Tel : 02-3461-XXXX


# Medical Conferences and PR Activities

Czy technologia Braster ma sens?  
Joanna Rusecka, 22.03.2019. Tagi: rak piersi, onkologia, termografia kontaktowa, termografia, braster

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